



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medicare Gynaecology Couche (SM8593D) (New)

- New,
- Made in UK,
- Electric control of height and tilt angle,
- Supplied with matching operators chair and head cushion,
- Compact frame design with retractable wheel system (RWD),
- Crescent shaped seat design,
- Basics:
 - Paper Roll holder (6005),
 - Removable headrest - pillow (6007),
 - Side support loops (6041),
 - Debris tray,
 - Removable patient leg support (6043),
 - Short frame cover (6052),
 - operators chair (6100),
- Base frame powder coated with antibacterial layer,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 178 cm,
 - Width: 65 cm,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from horizontal to 80°,
 - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
 - Fully adjustable patient leg supports (6043),
 - Battery back-up (6057),
 - Triple foot switch (6209),
- Technical Passport,
- 24-month warranty for electric



Fold-down leg extension



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD +[ZA DOPŁATK.]

