



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Shower Trolley (SM7700) (New)

- New
- Made in UK
- Large mobile wheel design with central locking and steering facility
- Hydraulic height range adjustment of 41cm (16") to 94cm (37")
- Foam padded PVC liner with deep sides and end panels for comfort and security
- Supplied with foam padded PVC patient support cushion
- The side support rails fold down to assist patient transfer
- Low gradient tilt towards foot end for water drainage via integral waste outlet
- Plug and drain hose supplied complete with storage bracket for hose
- Bumpers fitted to all corners for protection during transport procedures
- Bumpers can be removed to reduce the width to 65cm to access tight doorways
- Lowers to wheelchair height for ease of patient transfer
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Technical Data:
 - Length: 195 cm,
 - Width: 74 cm,
 - Safe working load and lifting capacity (from minimum height) of 240kgs (528lbs=37 stone)
 - Adjustable height: 41 cm to 94 cm,
- **Additional Equipment (paid):**
 - Spare / replacement PVC liner (6093)
 - Spare patient PVC cushion (6094)
- Technical Passport
- 24 months warranty for hydraulics

