



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medicare Bariatric 2 Section (SM4542) (New)

- New,
- Made in UK,
- Bariatric 2 section couch with retractable wheels,
- Remote Control,
- Vertical lifting design with stability throughout,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 200 cm,
 - Width: 105 cm,
 - Safe working load and lifting capacity of 325 kg (715 lbs),
 - Adjustable angle electric backrest from 0° to 75°,
 - Excellent height range from 45 cm to 98 cm,
- **Additional Equipment (extra paid):**
 - Breathing hole and plug for prone treatments (6004),
 - Paper Roll handle (6005),
 - Side support rails (Pair) (6074),
 - Side support rails cushions (Pair) (6074P),
 - Memory foam upholstery (60 mm) for additional patient comfort (6076),
 - Double Foot switch (6207),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD +[ZA DOPŁATKĄ]

