




Medicare Bariatric 3 Section (SM3573) (New)

- New,
- Made in UK,
- Bariatric 2 section couch with retractable wheels,
- Remote Control,
- Vertical lifting design with stability throughout,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 191 cm,
 - Width: 80 cm,
 - Safe working load and lifting capacity of 325 kg (715 lbs),
 - Adjustable angle electric backrest from 0° to 75 °,
 - Adjustable angle of footrest from -25° to -45°,
 - Excellent height range from 45 cm to 98 cm,
- **Additional Equipment (extra paid):**
 - Breathing hole and plug for prone treatments (6004),
 - Removable headrest (6007),
 - Paper Roll handle (6005),
 - Side support rails (Pair) (6074),
 - Side support rails cushions (Pair) (6074P),
 - Memory foam upholstery (60 mm) for additional patient comfort (6076),
 - Triple Foot switch (6209),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD (+ZA DOPŁAT)

