



Medicare Bariatric Drop End (SM3584) (New)

- New,
- Made in UK,
- Electric height elevation,
- 2 section couch with retractable wheels,
- Memory foam upholstery (60 mm) for additional patient comfort,
- Vertical lifting design with stability throughout,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Remote Control,
- Technical Specifications:

- Length: 180 cm,
- Width: 85 cm,
- Safe working load and lifting capacity of 325 kg (715 lbs),
- Adjustable angle electric backrest from 0° to 80°,
- Adjustable angle of footrest from 0° to -85°,
- Excellent height range from 50 cm to 105 cm,

• Additional Equipment (extra paid):

- Breathing hole and plug for prone treatments (6004),
- Removable headrest (6007)
- Paper Roll handle (6005),
- Side support rails (Pair) (6074),
- Side support rails cushions (Pair) (6074P),
- Technical Passport,
- 24-month warranty for electric

