



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl




Medicare Bariatric Drop End (SM3583) (New)

- New,
- Made in UK,
- Electric height, backrest and dual foot sections,
- 2 section couch with retractable wheels,
- Memory foam upholstery (60 mm) for additional patient comfort,
- Vertical lifting design with stability throughout,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Remote Control,
- Technical Specifications:

- Length: 180 cm,
- Width: 85 cm,
- Safe working load and lifting capacity of 325 kg (715 lbs),
- Adjustable angle electric backrest from 0° to 80°,
- Adjustable angle of footrest from 0° to -85°,
- Excellent height range from 50 cm to 105 cm,

• **Additional Equipment (extra paid):**

- Breathing hole and plug for prone treatments (6004),
- Removable headrest (6007),
- Paper Roll handle (6005),
- Side support rails (Pair) (6074),
- Side support rails cushions (Pair) (6074P),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPLATĄ)

