




Innovation Deluxe Dialysis (MG3490) (New)

- New,
- Made in UK,
- Electric Height Elevation with Electric Back & footrest and 2-way tilt - Individual locking castors,
- 2-way tilt ensures smooth transition when positioning patient for specialist procedures and treatments,
- Basic equipment:
 - Removable headrest (6007)
 - Paper Roll handle (6005),
- Compact frame design,
- Anti-bacterial powder coated framework,
- Acrylic capped vac formings to facilitate easy cleaning SEERS ComfortPlus™,
- Vertical lifting design with stability throughout,
- Fitted with multi adjustable armrests
- Fitted with adjustable patient footboard
- Fitted with pushing handles
- TiMotion electric linear actuators,
- Adjustable angle:
 - Trendelenburg position,
 - Auto-electric CPR functionality and auto-levelling software returns the couch quickly to a fully horizontal position ,
- Technical Specifications:
 - Length: 199 cm,
 - Width: 70 cm made from 60 mm foam,
 - Safe working load and lifting capacity of 250 kg (550lbs),
 - Adjustable angle electric backrest from 0 ° to 80 °,
 - Adjustable angle of footrest from - 35 ° to 10 °,
 - Excellent height range from 47 cm to 97 cm,
 - Wheels Tente® Ø125mm,
- **Additional Equipment (paid):**
 - Battery back-up (6057),
 - IV Pole (6054),
 - Tray (6241)
- Technical Passport,
- 24 months warranty for electric



TAPICERKI STANDARD [W CENIE]

