



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl


Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Phlebotomy Couch (SM9586P) (New)

- New
- Made in UK
- Electric with electric backrest and tilt
- Patient handgrip supports – to aid mount/dismount from the couch and provide patient security
- Fully adjustable phlebotomy arms for use on either side of couch
- Fitted with paper roll holder
- Compact frame design with retractable wheel system
- Adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- The seat section is inclined for patient security
- Clearance beneath frame for mobile hoist
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 180 cm,
 - Standard width 65cm (25")
 - Safe working load and lifting capacity (from minimum height) of 240kgs (528lbs = 37.7 stone)
 - Adjustable angle backrest from -10 to +85 degrees positive
 - Adjustable angle feet plates 0 ° to -85 °,
 - Excellent height range from 50-102cm (19"-39")
- **Additional Equipment (paid):**
 - Removable head cushion (6007)
 - Protective base frame shrouding (6052)
 - Battery back-up (6057),
 - High Operators Chair (6101)
 - Deluxe Saddle Seat (6107)
 - Standard Operators Stool (6103),
 - Triple foot switch operation (6200)
- Technical Passport
- 24 months warranty for electrics and gas springs



TAPICERKI STANDARD (W CENIE)

