



Examination Table Therapy 2 Section (ST2567) (New)

- New,
- Made in UK,
- Electric height elevation,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Hand switch fitted as standard on electric models,
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 189 cm,
 - Width: 63 cm,
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
 - Adjustable headrest: -25° to 75°,
 - Excellent height range from 45 cm (18") to 98 cm (38"),
- **Additional Equipment (extra paid):**
 - Perimeter foot switch (EPFS) (ST6001),
 - Side Support Rails (ST6002),
 - Side Support Cushions (ST6002P),
 - Paper towel holder (ST6005),
 - Upper limb support (ST6006),
 - Removable head/knee support (ST6007),
 - Matching therapy chair or stool
 - Single foot switch (ST6200),
 - Double foot switch (ST6239),
- Technical Passport,
- 24-month warranty for electrics and gas springs

