



## Examination Table Therapy 3 Section (ST3559) (New)

- New,
- Made in UK,
- Hydraulic height elevation,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Plus head design (with adjustable height armrests),
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- TiMotion electric linear actuators,
- Technical Data:
  - Length: 188 cm,
  - Width: 65 cm,
  - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
  - Adjustable headrest: -25° to 75°,
  - Adjustable Backrest: 0° do 80°,
  - Excellent height range from 45 cm (18") to 98 cm (38"),
- **Additional Equipment (extra paid):**
  - Side Support Rails (ST6002),
  - Side Support Cushions (ST6002P),
  - Paper towel holder (ST6005),
  - Upper limb support (ST6006),
  - Removable head/knee support (ST6007),
  - Matching therapy chair or stool,
- Technical Passport,
- 24-month warranty for electrics and gas springs

