



Examination Table Deluxe Therapy Drainage (ST3547S) (New)

- New,
- Made in UK,
- Electric height elevation via integral perimeter foot switch,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- Deluxe Therapy Couch with Drainage & Plus head design - split leg,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 208 cm,
 - Width: 65 cm,
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
 - Adjustable headrest: -25° to 40°,
 - Adjustable leg plates: 0° to 90°,
 - Excellent height range from 46 - 99 cm (19"-39")
- **Additional Equipment (extra paid):**
 - Paper towel holder (ST6005),
 - Removable head/knee support (ST6007),
 - Protective Base Shroud (ST6053),
 - Matching therapy chair or stool,
- Technical Passport,
- 24-month warranty for electrics and gas springs

