



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl


Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Table Deluxe Therapy Non - Drainage (ST3349) (New)

- New,
- Made in UK,
- Electric height elevation via integral perimeter foot switch,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- Deluxe Therapy Couch with head design, Perimeter foot switch operation for electric models (ST6001) - armrests,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 195 cm,
 - Width: 65 cm,
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
 - Adjustable headrest: -25° to 40°,
 - Adjustable leg plates: 0° to 90°,
 - Excellent height range from 46 - 99 cm (19"-39"),
- **Additional Equipment (extra paid):**
 - Paper towel holder (ST6005),
 - Removable head/knee support (ST6007),
 - Protective Base Shroud (ST6053),
 - Matching therapy chair or stool,
- Technical Passport,
- 24-month warranty for electrics and gas springs



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD +[ZA DOPŁATĄ]

