



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl


Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Table Therapy Traction Couches and Packages (ST6557) (New)

- New,
- Made in UK,
- Hydraulic Height Elevation,
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Head section fitted with breathing hole and plug,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 188 cm,
 - Width: 70 cm,
 - Safe working load and lifting capacity (from minimum height) of 240 kg (528 lbs = 37 stone),
 - Positive & Negative head design with adjustable angle of -25 to +75 degrees,
 - Adjustable Backrest: 0° to 80°,
 - Excellent height range from 50 - 101 cm (20" - 40"),
- **Additional Equipment (extra paid):**
 - Pelvic (Lumbar) traction harness (ST6046),
 - Thoracic traction harness (ST6047),
 - Thoracic harness straps (pair) (ST6048),
 - Cervical traction head halter (ST6049),
 - Cervical pillow (ST6067),
 - Flexion Stool (ST6069),
 - Spreader bar (ST6070),
 - Arm for neck accessories (ST6071),
 - Traction SM12 (ST6500),
- Technical Passport,
- 24-month warranty for electrics and gas springs



TAPICERKI STANDARD [W CENIE]

