



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Table Therapy Bobath Couches (ST4542W) (New)

- New,
- Made in UK,
- Electric height elevation,
- Retractable wheel design and adjustable foot for stability on uneven surfaces,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Hand switch fitted as standard on electric models,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 200 cm,
 - Width: 125 cm,
 - Safe working load and lifting capacity (from minimum height) of 250 kg (550 lbs = 39 stone),
 - Adjustable angle backrest facility (horizontal to 75 degrees),
 - Excellent height range from 46 cm (18") to 100 cm (39"),
- **Additional Equipment (extra paid):**
 - Perimeter foot switch operation for electric models (EPFS) (ST6001),
 - Breathing hole and plug for prone treatments (ST6004),
 - Side Support Rails (each) (ST6074),
 - Side Support Rail Cushion (each) (ST6074P),
 - Single foot switch (ST6200),
 - Double foot switch (ST6239),
- Technical Passport,
- 24-month warranty for electrics and gas springs

