




Examination Table Therapy Deluxe Bobath Couches (ST4642) (New)

- New,
- Made in UK,
- Electric height operation,
- Large wheel design with steering facility and central locking from a single level operation,
- Hand switch fitted as standard,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 200 cm,
 - Width: 105 cm,
 - Safe working load and lifting capacity (from minimum height) of 325 kg (51 stone),
 - Powered electric backrest fully adjustable from horizontal to 70 degrees positive,
 - Excellent height range of 45 - 98 cm (18-38"),
- **Additional Equipment (extra paid):**
 - Bobath patient limb support for use only at food end of couch (ST6012),
 - Patient handgrip supports for use only at food end of couch (pair) (ST6013),
 - Accessory rail for fixation on the couch - pair (ST6014),
 - Battery back-up facility (ST6057),
 - Accessory rail for fixation of patient handgrip supports on the side of couch - one required per side (ST6073),
 - Side Support Rails (each) (ST6074),
 - Side Support Rail Cushion (each) (ST6074P),
 - Single foot switch (ST6200),
 - Double foot switch (ST6239),
- Technical Passport,
- 24-month warranty for electrics and gas springs



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATĄ)

