



Examination Table Therapy Tilt Table (ST7641DL) (New)

- New,
- Made in UK,
- Electric height and tilting operation with hand switch control,
- Large wheel design with central locking and steering facility on Standard model,
- Individual braking castors on Economy model,
- Lowers to wheelchair height for ease of patient transfer,
- Clearance beneath frame for mobile hoist,
- Remote Control,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 195 cm,
 - Width: 70 cm,
 - Safe working load and lifting capacity (from minimum height) of 225 kg (495 lbs),
 - Height range 47 cm to 101 cm (18 - 39"),
 - Tilt angle from horizontal to 90 degrees vertical,
- **Additional Equipment (extra paid):**
 - Removable headrest (ST6007),
 - Battery back-up facility (ST6057),
 - Hip/Leg Stabilisers (pair) (ST6058),
 - Adjustable Pommel Support (ST6059),
 - Child Footboard Facility (ST6060),
 - Dual Leg Harness Set (ST6063),
 - Pair Dual Knee Harnesses (ST6064),
- Technical Passport,
- 24-month warranty for electrics and gas springs

