



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl


Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rehabilitation table Therapy Hygiene Table (ST1551-SEERSMEDICAL) (New)

New,
English production,
Table height adjusted by means of a hydraulic cylinder,
Retractable wheels and adjustable feet provide excellent stability (RWD),
Frame construction designed for easy use of a mobile lift,
Cylinders used in the TiMotion table,
Technical data:
Total length of the table bed: 186 cm,
The width of the bed is 65 cm,
Safe working load and lifting capacity (from the minimum height) 240 kg (528 lbs)
Table height adjustment range from 46 cm to 98 cm,
Additional accessories (extra charge):
Rail side barriers - pair (ST6002),
Upholstered side barriers - pair (ST6002P),
Paper hanger (ST6005),
Cover with a breathing hole (ST6010),
Cover without a breathing hole (ST6011),
Removable headrest - Cushion (ST6007),
Table base protective cover (ST6053),
Matching chair or driver's stool,
Has a current Technical Passport issued,
Warranty 24 - months for hydraulic cylinders



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATĄ)

