



Examination Table Therapy Hygiene Table (ST1651B) (New)

- New,
- Made in UK,
- Hygiene Therapy Table - hydraulic, large mobile wheel design,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- TiMotion electric linear actuators,
- Technical Data:
 - Length: 152 cm,
 - Width: 65 cm,
 - Safe working load and lifting capacity (from minimum height) of 240 kgs (528 lbs = 37 stone),
 - Excellent height range from 43 cm (17") to 97 cm (37"),
- **Additional Equipment (extra paid):**
 - Side support rails (ST6002),
 - Side support cushions (ST6002P),
 - Paper towel holder (ST6005),
 - Stretch towelling covers (ST6010), (ST6011),
 - Removable head/knee support (ST6007),
 - Protective base shrouding (ST6053),
 - Matching therapy chair or stool,
- Technical Passport,
- 24-month warranty for electrics and gas springs

