



Examination Table Sterling 3 Section Couch (SX3008) (New)

- New,
- Made in UK,
- Electric Elevation,
- 100 mm diameter individual locking castors,
- Remote Control,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 188 cm,
 - Width: 63 cm,
 - Safe working load and lifting capacity (from minimum height) of 200 kgs (440 lbs = 31.5 stone),
 - Adjustable angle electric backrest from 0° to 85°,
 - Adjustable angle of footrest/headrest from -20° to 40°,
 - Excellent height range from 46 cm to 99 cm,
- **Additional Equipment (extra paid):**
 - Side support rails (Pair) (6003),
 - Side support rails cushions (Pair) (6003P),
 - Breathing hole and plug for prone treatments (6004),
 - Removable headrest (6007),
 - Paper Roll handle (6005),
 - Foot switch (6200),
- Technical Passport,
- 24-month warranty for electrics and gas springs

