



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Saddle Stool (MC6145) (New)

- New
- Made in UK
- The curved seat pad design provides a similar ergonomic seated position to the general saddle stool
- All models feature a polished aluminium five star base and castors.
- Saddle stool - high model
- Height range 54cm to 74cm
- Foot support ring fitted as standard.
- Max lifting capacity: 150kg,
- 24 months warranty



TAPICERKI STANDARD + (ZA DOPŁATĄ)



TAPICERKI STANDARD



TAPICERKI STANDARD

TAPICERKI STANDARD (W CENIE)

