



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Contoured Medical Chair (MC7212) (New)

- New
- Made in UK
- The 3 lever mechanism provides independent positioning of the backrest angle, seat pad tilt and chair height
- Contoured medical chair - standard model
- Height range 45cm to 59cm
- Max lifting capacity: 150kg,
- 24 months warranty



TAPICERKI STANDARD + [ZA DOPŁATK.]



TAPICERKI STANDARD (W CENIE)

