



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Contoured Medical Chair (MC7213) (New)**

- New
- Made in UK
- A versatile chair with a contoured seat pad and raised front to promote positive posture whilst seated
- Contoured medical chair - high model
- Height range 54cm to 74cm
- Foot support ring fitted as standard.
- Max lifting capacity: 150kg,
- 24 months warranty



TAPICERKI STANDARD + [ZA DOPŁATK.]



TAPICERKI STANDARD (W CENIE )

