



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Care Chair General Medical (MC6178) (New)

- New
- Made in UK
- Ideal for all bedside requirements in NHS and nursing home environments
- Various waterproof and anti-bacterial vinyl and fabric options available
- Removable and reversible seat cushions
- Easy to clean and wipe-able
- Anti-bacterial lacquered wooden frames
- Strong and durable wooden design
- Max lifting capacity: 150kg,
- Adjustable height: 45 - 59 cm,
- **Additional Equipment (paid):**
 - A set of anti-static castors (6178),
 - Load loading castors automatically locked (6455),
 - Manual Locking castors (6172),
 - Drop down arm supports (6175),
 - Foot Height Control (FHC) (6174),
 - Foot support ring (6176),
- 24 months warranty



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATĄ)

