



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Chair Surgeons & Sonographers (MC6158) (New)

- New
- Made in UK
- Surgeons chair - standard model
- Height range 45cm to 59cm.
- Wide choice of models available
- Ergonomically designed seat
- Compact backrest for improved scanning procedures
- Seat and backrest upholstered on both sides
- Weight capacity and safe working load of 150Kg
- Anti-bacterial and waterproof vinyl upholstery
- Restraint castors fitted as standard
- **Dostępne akcesoria (za dopłatą):**
 - A set of anti-static castors (6178),
 - Load loading castors automatically locked (6455),
 - Manual Locking castors (6172),
 - Glides (static feet) (6173),
 - Drop down arm supports (6175),
 - Foot Height Control (FHC) (6174),
 - Foot support ring (6176),
- 24 months warranty



TAPICERKI STANDARD + (ZA DOPŁATĄ)



TAPICERKI STANDARD (W CENIE)

