



Clinical Chair Saddle (MC6113) (New)

- New
- Made in UK
- Clinical Chair with Arms
- Five star base and castors
- Ergonomic seat shape
- Max lifting capacity: 150kg,
- Adjustable height: 41 - 55 cm,
- **Additional Equipment (paid):**
 - A set of anti-static castors (6178),
 - Load loading castors automatically locked (6455),
 - Manual Locking castors (6172),
 - Drop down arm supports (6175),
 - Foot Height Control (FHC) (6174),
 - Foot support ring (6176),
- 24 months warranty

