



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Clinical Chair (MC 6115) (New)

- New
- Made in UK
- Clinical Chair with Arms
- Five star base and castors and ergonomic seat shape
- Max lifting capacity: 150kg,
- Adjustable height: 50 – 70 cm,
- Foot support ring as standard
- **Additional Equipment (paid):**
 - A set of anti-static castors (6178),
 - Load loading castors automatically locked (6455),
 - Manual Locking castors (6172),
 - Drop down arm supports (6175),
 - Foot Height Control (FHC) (6174),
- 24 months warranty



MC6115
54cm to 74cm



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATK.)

