



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Multi Procedures Chair (MC6163) (New)

- New
- Made in UK,
- Chair with 360° swing around arm - High model
- Height range 54cm to 74cm
- Foot support ring fitted as standard.
- Max lifting capacity: 150kg,
- 24 months warranty



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATĄ)

