



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Examination Chair SEERSMEDICAL (MC6172) (New)

- New
- Made in UK,
- Height regulated by gas spring,
- Removable headrest,
- Pentacle frame,
- Maximum load: 150 kg,
- 24-month warranty



### TAPICERKI STANDARD (W CENIE )



### TAPICERKI STANDARD + (ZA DOPŁATK.)

