



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Phlebotomy Chair (MC6154) (New)

- New
- Made in UK
- The chair is a compact and space saving alternative to electrically operated chairs and couches.
- Supplied with a sturdy base on glides (static feet) or with several castor options.
- Integral arm supports are fitted as standard to assist the patient in mounting and dismounting the chair.
- The chair has a height range of 49cm to 63cm
- Max load: 200kg,
- 24 months warranty



TAPICERKI STANDARD + [ZA DOPŁATK.]



TAPICERKI STANDARD (W CENIE)

