



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient Chair 5106 (New)

- New,
- Made in Germany,
- Universal patient chair suitable for ENT and other medical examinations,
- With infinite height adjustment using a double pedal switch,
- Maintenance-free, electromotive lifting equipment, lift height 200 mm with 12 V-control in accordance with IEC 601-1,
- The seat is horizontally mounted,
- The back rest is adjustable from the vertical (95°) to below the horizontal position by means of a gas spring,
- The headrest can be infinitely adjusted in height and, if necessary, can be pulled all the way out,
- The two armrests can be folded down towards the back and the front,
- The entire upper part can be fully rotated and locked in any position,
- The one-piece upholstery is secured with Velcro straps,
- **Additional Equipment (extra paid):**
 - Foot switch,
 - Headrest with dual connection,
 - Headrest,
 - ENT,
- Weight: 80 kg,
- Technical Passport,
- 24-month warranty

