



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Patient Chair JORG & SOHN 5101/3 (New)

- New
- Made in Germany,
- Reliable patient chair for general use, particularly suitable for ENT and ophthalmology departments,
- Large base plate,
- Maintenance-free electro-mechanical lifting equipment, lift height 300 mm,
- Seat height: 51 - 81 cm,
- Height adjustment with flat foot pedal,
- Comfortably padded top can be rotated by 360° and lockable in any position,
- The high back rest is adjustable from vertical to horizontal position by 90° by means of a gas spring,
- Pair of height adjustable padded armrests which can be lowered below the seat cushion height (so the patient can be placed into the chair at 90°),
- Foot rest with fold out footboard can be tilted by 90° and rotated to the upward horizontal position, can be pulled out by 18 cm. The chair is convertible into a sturdy bed,
- Removable and inclinable headrest,
- Upholstery as per our sample card,
- Painted metal components RAL 9002,
- **Additional Equipment (extra paid):**
 - Neutral position - pedal switch,
 - Headrest with dual switch,
 - Headrest,
 - Footrest,
 - Regulated armrests,
- Technical Passport,
- 24-month warranty

