



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Patient Chair JORG & SOHN 5104/3 (New)**

- New
- Made in Germany,
- Universal patient chair suitable for ENT and other medical examinations,
- Height is infinitely adjustable using a double pedal switch,
- Maintenance-free, electromotive lifting equipment, lift height 200 mm with 12 V-control,
- The seat is horizontally mounted (slight seat inclination),
- The back rest is adjustable from the vertical (95°) to below the horizontal position by means of a gas spring,
- The headrest can be infinitely adjusted in height and, if necessary, can be pulled all the way out. After loosening the hand screw, it can be additionally adjusted by approx. 170° in inclination,
- The entire leg rest can be folded into an upward horizontal position. The leg rest is extensible by synchronously revolving footrest,
- The two comfortably upholstered armrests can be folded down towards the back and the front,
- The entire upper part can be fully rotated and locked in any position,
- The one-piece upholstery is secured with Velcro straps,
- Weight: 80 kg,
- **Additional Equipment (extra paid):**
  - Neutral position - pedal switch,
  - Headrest with dual switch,
  - Headrest,
  - Footrest,
  - ENT,
- Technical Passport,
- 24-month warranty

