



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Patient Chair JORG & SOHN 5106 (New)**

- New
- Made in Germany,
- Universal patient chair suitable for ENT and other medical examinations,
- With infinite height adjustment using a double pedal switch,
- Maintenance-free, electromotive lifting equipment, lift height 200 mm with 12 V-control in accordance with IEC 601-1,
- The seat is horizontally mounted,
- The back rest is adjustable from the vertical (95°) to below the horizontal position by means of a gas spring,
- The headrest can be infinitely adjusted in height and, if necessary, can be pulled all the way out,
- The two armrests can be folded down towards the back and the front,
- The entire upper part can be fully rotated and locked in any position,
- The one-piece upholstery is secured with Velcro straps,
- **Additional Equipment (extra paid):**
  - Foot switch,
  - Headrest with dual connection,
  - Headrest,
  - ENT,
- Weight: 80 kg,
- Technical Passport,
- 24-month warranty

