



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Chair Coburg Ray-O-Seat JORG & SOHN 4045 U (New)

- New,
- Made in Germany,
- X-ray, imprint and multi-purpose chair suitable for orthodontics and prophylaxis with hydraulic pump lever,
- Large cast aluminium base plate 50 x 50 cm,
- Standard version with 4 fixed rollers (move in one direction without brake),
- Available with feet instead of rollers if required,
- Anthracite plastic bucket seat,
- Metal parts powder-coated silver and black,
- Adjustments options:
 - lift height 19 cm,
 - seat height 40 - 59 cm,
 - upper part can be rotated by 360° and locked,
 - optional headrest fully adjustable,
- **Additional Equipment (extra paid):**
 - cushion pad for seat (as per our sample cards),
 - ring armrests,
 - shell headrest, black polyurethane,
 - shell headrest, upholstered,
- Technical Passport,
- 24-month warranty

