



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Chair Coburg Ray-O-Seat JORG & SOHN 4046 EH-U (New)

- New,
- Made in Germany,
- X-ray, impression and multipurpose chair for orthodontics and prophylaxis, with electromechanical lift,
- Height adjustment by foot pedal,
- Large base plate made of aluminium cast 50 x 50 cm,
- 4 small casters (plastic studs on request),
- Anthracite plastic bucket seat,
- Metal parts painted in RAL 9002 and RAL 9005,
- Adjustments options:
 - stroke: 20 cm,
 - height of seat: 51 - 71 cm,
 - optional headrest fully adjustable,
- **Additional Equipment (extra paid):**
 - cushion pad for seat (as per our sample card),
 - rotating device (upper part can be rotated by 360° and locked),
 - ring armrests,
 - shell headrest, black polyurethane,
 - shell headrest, upholstered (as per our sample card),
 - Footstool,
- Technical Passport,
- 24-month warranty

