



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Chair Coburg X-Ray-Lift JORG & SOHN 4047 (New)

- New,
- Made in Germany,
- X-ray chair with 5-beam cast aluminium safety base,
- Upper part can be rotated,
- Height is adjusted by hand-operated safety gas spring,
- With 5 special rollers which lock once the chair is loaded,
- Modern anthracite plastic seat frame,
- Powder-coated silver and black metal components,
- Adjustment Options:
 - lift height 20 cm,
 - seat height 45 - 65 cm,
 - optional headrest fully adjustable,
- **Additional Equipment (extra paid):**
 - upper part cannot be rotated,
 - cushion pad for seat (as per our sample card),
 - ring armrests,
 - shell headrest, black polyurethane,
 - shell headrest, upholstered,
- Technical Passport,
- 24-month warranty

