



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **X-ray Wall-mounted Chair Coburg RoWa 4055 (New)**

- New
- Made in Germany
- elegant, wall-mounted X-ray chair, a further development of the approved RÖWA 4050
- chair with a stable base support
- anatomically padded seat and back rest
- plastic cover on seat base
- seat can be folded away and is automatically locked in sitting or reclining position
- integrated profile rail in the back rest enables precise positioning of the headrest
- upholstery as per our sample card
- powder-coated metal components RAL 9002
- **Additional equipment (paid):**
  - wedge-shaped headrest
  - u-shaped headrest with self-locking slider for the profile rail (single-handed operation)
  - height adjustable radiation protection shield with wall mounting attachment
  - seat cushion for children
- Technical Passport
- 24 months warranty

