



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



X-ray Wall-mounted Chair Coburg RoWa 4050 (New)

- New
- Made in Germany
- world-renowned, wall-mounted X-ray chair, space-saving and practical
- seat with stabile floor-mounting support
- tiltable seat
- seat and backrest with ergonomic upholstery as per our sample card
- backrest with integrated Velcro fastener for fully adjustable headrest
- metal parts powder-coated RAL 9002
- **Additional equipment (paid):**
 - wedge-shaped headrest
 - U-shaped headrest
 - height adjustable X-ray protection shield with wall mounting attachment
 - seat cushion for children
- Technical Passport
- 24 months warranty

