



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Dentalift 2515 (New)

- New
- Made in Germany
- dental operator's stool (according to DIN EN ISO 7493)
- chrome safety base with 5 twin castors
- chrome foot ring
- anatomically-shaped seat- and back rest
- special seat mechanism for individual adjustment of seat and back rest
- upholstery colours as per our sample cards
- Adjustment options:
 - moulded seat infinitely adjustable in inclination
 - back rest tilt mechanism permanently oscillation and can be locked into position
 - back rest height can be adjusted by 7.5 cm
 - infinite seat height adjustment by means of a safety gas spring and foot control disc
- Seat Height:
 - L: 43 - 57 cm,
 - T: 49 - 69 cm,
- **Additional Equipment (paid):**
 - upholstery colours as per our sample cards
 - left arm
 - right arm
- Technical Passport
- 24 months warranty

