



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67  
 E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



**Medical Stool Coburg Dentalift 22015 (New)**

- New
- Made in Germany
- dental operator's stool (according to DIN EN ISO 7493)
- base made of stainless steel with 5 twin castors for hard floors
- chrome-plated foot ring
- height adjustment using a safety gas spring by means of an overlying foot-controlled disc
- seat and back rest anatomically-shaped and easily adjustable by means of a lever
- comfortable upholstery as per our sample cards (special colours subject to a surcharge)
- Seat Height:
  - L: 42 - 56 cm,
  - T: 48 - 68 cm,
- **Additional Equipment (paid):**
  - metal parts in WHITE LINE - RAL 9002, powder-coated
  - operator's armrest 2014, swivel and height adjustable, mountable at right or/and left side
  - special heights on request
- Technical Passport
- 24 months warranty

