



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Dentalift 11001 (New)

- New
- Made in Germany
- dental operator's stool (according to DIN EN ISO 7493)
- safety base made of stainless steel with 5 twin castors for hard floors
- chrome-plated foot ring height adjustment using a safety gas spring controlled by means of a foot-controlled disc comfortable upholstery as per our sample cards (special colours subject to a surcharge)
- Airflow rounded seat Ø38cm
- Seat Height:
 - low version: 44 - 57 cm,
 - high version: 51 - 70 cm,
- **Additional Equipment (paid):**
 - back rests no. 1 / 2 / 3 / 4
 - moulded seat 1103
 - Super-Denta
 - Airflow rounded seat Ø35cm
 - operator's armrest 2014, swivel and height adjustable, mountable at right or/and left side
 - metal parts in WHITE LINE (RAL 9002), powder-coated
 - seat angle adjustment with threaded spindle for individual adjustments of the seat inclination in moulded seats
- Technical Passport
- 24 months warranty

