



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Dentalift 11004 (New)

- New
- Made in Germany
- dental operator's stool (according to DIN EN ISO 7493)
- safety base made of stainless steel with 5 twin castors for hard floors
- anatomically-shaped padded seat infinite height adjustment using a hand lever and a safety gas spring
- seat angle adjustment with threaded spindle for individual adjustment of the seat inclination
- upholstery colours as per our sample cards
- special colours subject to a surcharge
- seat height: 44 - 57 cm
- **Additional Equipment (paid):**
 - operator's armrests 2014, swivel and height adjustable, mountable at right or/and left side
 - special colours subject to a surcharge
 - metal parts in WHITE LINE - RAL 9002, powder-coated
 - castors for soft floors
- Technical Passport
- 24 months warranty

