



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Dentalift 11005 (New)

- New
- Made in Germany
- dental operator's stool (according to DIN EN ISO 7493)
- safety base made of stainless steel with 5 twin castors (secured by extended bolt) for hard floors
- infinite height adjustment using a hand lever and a safety gas spring
- only available with Airflow rounded seat Ø35cm and back rest no. 3
- comfortable upholstery as per our sample cards /special colours subject to a surcharge
- seat height: 56 - 75cm
- **Additional Equipment (paid):**
 - special colours subject to a surcharge
 - metal parts in WHITE LINE - RAL 9002, powder-coated
 - castors for soft floors
- Technical Passport
- 24 months warranty

