



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Dentalift 22006 (New)

- New
- Made in Germany
- dental operator's stool (according to DIN EN ISO 7493)
- safety base made of stainless steel with 5 twin castors for hard floors
- Airflow rounded seat Ø35cm
- infinite height adjustment by means of a hand lever and a safety gas spring
- comfortable upholstery as per our sample cards
- special colours subject to a surcharge
- Seat height:
 - L: 42 - 55 cm,
 - T: 46 - 65 cm,
- **Additional Equipment (paid):**
 - special colours subject to a surcharge
 - back rests no. 1 / 2 / 3 / 4
 - moulded seat 1103
 - Super-Denta moulded seat
 - Airflow round seat Ø38cm
 - metal parts in WHITE LINE - RAL 9002, powder-coated
 - operator's armrests 2014, swivel and height adjustable, mountable at right or/and left side
 - seat angle adjustment with threaded spindle for the individual adjustment of the seat inclination in moulded seats
 - castors for soft floors
- Technical Passport
- 24 months warranty

