



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybierański
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Eurolift 11015 (New)

- New
- Made in Germany
- multi-purpose stool (in accordance with DIN EN ISO 7493)
- safety base made of chrome-plated section steel tube with 5 twin rollers for hard floors
- chrome-plated foot ring
- height adjustment using a safety gas spring by means of a hand lever
- seat and back rest anatomically-shaped and easily adjustable by means of a lever
- comfortable upholstery as per our sample card
- special colours subject to a surcharge
- Seat height: 47 - 66 cm,
- **Additional Equipment (paid):**
 - special colours subject to a surcharge
 - adjustable armrests: left side/ right side mounted
- Technical Passport
- 24 months warranty

