



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Eurolift 33010 (New)

- New
- Made in Germany
- multi-purpose stool
- safety base made of chrome-plated steel with 5 twin rollers for hard floors
- seat and back rest anatomically-shaped and asynchronously infinitely adjustable
- infinite seat height adjustment using a hand lever and a safety gas spring
- upholstery colours as per our sample card
- special colours subject to a surcharge
- Seat Height:
 - L: 41 - 54 cm,
 - T: 46 - 65 cm,
- **Additional Equipment (paid):**
 - special colours subject to a surcharge
 - Arm support
- Technical Passport
- 24 months warranty

