



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Laborlift 532 (New)

- New
- Made in Germany
- quality laboratory stool
- large aluminium safety base with 5 twin castors
- infinite seat height adjustment by means of a hand lever and a safety gas spring
- back rest easily adjustable in height and depth using knob:
 - height = 7 cm
 - depth = 5.5 cm
- seat and back rest anatomically-padded
- comfortable upholstery as per our sample cards
- special colours subject to a surcharge
- seat height: 44 - 63 cm
- **Additional Equipment (paid):**
 - special colours subject to a surcharge
- Technical Passport
- 24 months warranty

