



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Ergolift 2520 (New)

- New
- Made in Germany
- new multi-purpose stool with ergonomic saddle seat
- stainless steel safety base with 5 twin castors for hard floors
- height adjustment using a safety gas spring controlled by means of a stainless foot control disc
- durable seat cover made from breathable synthetic leather as per our sample cards
- special colours subject to a surcharge
- seat height: 51 - 76cm
- **Additional Equipment (paid):**
 - special colours subject to a surcharge
 - Mini Seat (30 x 38 cm),
 - seat angle adjustment with threaded spindle for individual adjustment of seat inclination
- Technical Passport
- 24 months warranty

