



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Medicalift 3014 (New)

- New,
- Made in Germany,
- Special surgical stool suitable for all areas of surgery, particularly microsurgery (in accordance with DIN EN ISO 7493),
- Stainless steel safety base with 5 anti-static twin rollers,
- Chrome foot ring,
- Upper part can be rotated and locked into position,
- Seat and back rest anatomically shaped and easily adjustable by means of a lever,
- Comfortable upholstery, anti-static, black,
- Height adjustment via hydraulic pump lever; lift height: 19 cm,
- Seat Height: 53 - 72 cm,
- **Additional Equipment (extra paid):**
 - central stop bar,
 - operator armrests 2016,
 - 1 pair of flat stainless steel armrests, can be sterilized,
 - adjustable height and can be rotated on all sides by means of a ball joint,
 - available with washable upholstery if required,
- Technical Passport,
- 24-month warranty

