



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool Coburg Medicalift 3030 (New)

- New,
- Made in Germany,
- Special surgical chair suitable for all areas of surgery, particularly microsurgery,
- Stainless steel safety base with 5 anti-static twin rollers,
- Chrome foot ring,
- Upper part can be rotated and locked into position,
- Back rest can be adjusted by 7 cm,
- Comfortable upholstery, anti-static, black,
- Height adjustment with maintenance-free hydraulic pump lever, lifting height: 19 cm,
- Seat height: 53 - 72 cm,
- **Additional Equipment (extra paid):**
 - pair of upholstered operator armrests 2030,
 - central stop bar,
- Technical Passport,
- 24-month warranty

